

# AMBULANCEMED

AmbulanceMed, who combines technical expertise with the commercial ones and maintains with national/international experience, has led the delivery of thousands of ambulances until today.

## FACTORY - TURKEY

info@ambulancemed.com  
sales@ambulancemed.com

+90 554 969 50 08

DAĞYAKA MAH. 2018 SK. NO:23  
AKSAN TİCARET MERKEZİ  
KAHRAMANKAZAN / ANKARA

## FACTORY - DUBAI / UAE

info@ambulancemed.com  
dubai@ambulancemed.com

+971 56 332 6520

X3WR+Q57, Mina Jebel Ali, N506  
Jebel Ali Free Zone, Dubai, UAE

## OFFICE - DUBAI / UAE

info@ambulancemed.com  
dubai@ambulancemed.com

+90 554 969 50 08

Unit no. 604-605, 6F Opal Tower,  
Business Bay/ Dubai/ 128846  
Registration No: ICC20210433

## OFFICE / INTERNATIONAL

info@ambulancemed.com  
export@ambulancemed.com

+90 554 969 50 08

ATB İş Merkezi, G Blok No: 172  
Yenimahalle / Ankara



Ambulancemed

# CONTAINER MATERNITY UNIT

Rapid Response. Unmatched Care  
Container Clinics...







# *TABLE OF CONTENT*

You can explore a wealth of information about Container Maternity Unit within this brochure. Please refer to the following page numbers to find detailed content on each topic

01. [Container Maternity Unit](#)

02. [Vehicle Interior](#)

03. [Table Of Content](#)

04. [Technical Specifications](#)

05. [Medical Specifications](#)

06. [Gallery](#)

07. [Technical Equipments](#)

08. [Medical Equipments](#)

You are about to embark on an informative journey delving deep into the realm of Maternity Unit

# *WELCOME* INTRODUCTION

Immerse yourself in the realm of AmbulanceMed and discover the pinnacle of excellence in emergency healthcare solutions. Introducing our visionary company's latest achievement: the Container Maternity Unit – a masterpiece meticulously fashioned to elevate the benchmarks of medical mobility. Within this enchanting brochure, we extend an invitation to delve into the details of the AmbulanceMed creation, where aesthetics seamlessly harmonize with utility. Uncover the unparalleled sophistication of our state-of-the-art features, ingeniously designed to provide uncompromising maternal care and ensure rapid, effective emergency response. Embark with us on a journey that redefines innovation, propelling the evolution of future-ready maternity mobil solutions.



# CONTAINER MATERNITY UNIT

AmbulanceMed is dedicated to a transformative mission. Our newest innovation, the Container Maternity Unit, stands as a testament to this commitment. Crafted meticulously within our advanced Dubai and Ankara facilities, this unit is meticulously structured to offer unmatched medical excellence. Our design prioritizes impeccable hygiene, ensuring utmost comfort and care for both expectant mothers and medical staff.

Central to the Container Maternity Unit's innovation is its thoughtfully designed patient chamber, expertly tailored to uphold the most exacting medical norms. Recognizing the vital significance of maternity care, we have dedicated abundant space for optimal patient comfort and treatment. With meticulous precision, every facet of this compartment is purposefully conceived, accommodating seamless organization of specialized medical equipment, ensuring swift access during crucial moments in maternal emergencies.



## ■ *Advanced Illumination and Special Alert System*

The Container Maternity Unit introduces an energy-efficient lighting framework, ensuring optimal brightness for maternity procedures while prioritizing minimal power usage. Our specialized alert system is intricately tailored to capture attention while remaining energy-conscious. The harmonious integration of efficient lighting and alert mechanisms guarantees secure navigation through traffic, facilitating prompt arrivals during critical maternity situations.

At AmbulanceMed, our ethos embraces embedding cutting-edge technologies within our medical vehicles. The Container Maternity Unit have avant-garde technologies simplify the daily tasks of medical personnel, optimizing workflows, and elevating operational efficiency during maternity emergencies and procedures.

- The Container Maternity Unit, a creation of AmbulanceMed, transcends as a profound embodiment of our unwavering commitment to safety and well-being. It serves as a dynamic stage to unveil cutting-edge medical technologies, poised to redefine the horizon of maternal healthcare.



**AmbulanceMed**

As a premier innovator in medical mobility, our dedication transcends borders – reaching not only local communities but extending across the world. Our Container Maternity Unit, a testament to this commitment, finds its way to diverse nations, ensuring access to exceptional maternity care on a global scale.

The Container Maternity Unit exemplifies our resolute pursuit of excellence, safety, and forward thinking. With its unparalleled patient care amenities, ergonomic blueprint, efficient illumination and alert systems, it emerges as a beacon in the medical mobility landscape. Together, we strive to redefine maternity care standards, leaving an indelible mark on countless lives worldwide. Choose AmbulanceMed – your steadfast ally in delivering unparalleled medical mobility solutions.







## UNIT INTERIOR

AmbulanceMed embraces a profound sense of accomplishment in outfitting the interior of the Container Maternity Unit with an array of state-of-the-art technological marvels, in perfect harmony with the dynamic requirements of our contemporary world. Our unwavering commitment to elevating maternal healthcare to the highest echelons is impeccably manifested in every facet of this meticulously designed vehicle.

The interior space is meticulously appointed with cutting-edge, custom designed equipment, ensuring comprehensive readiness to handle diverse maternity scenarios. From avant-garde medical instruments tailored to maternity care to intuitively arranged layouts, each element has been meticulously curated to cater to the distinct demands of maternal healthcare practitioners and their patients. AmbulanceMed stands as the trusted vanguard, guaranteeing that our technologically advanced Container Maternity Unit will deliver unparalleled care for women's well-being, precisely when it is needed most.







## GENERAL DESCRIPTIONS

### CONTAINER MATERNITY CLINIC

Mobile Maternity clinic is formed by integrating ISUZU NPR truck and 20 ft container. It is equipped to be used for childbirth and gynecological examination.

The clinic contains, a gynecological examination chair, an examination bed and step, bassinet, storage areas, cabinets and drawers, outside area, bench tap, examination lamp, crash card, air conditioning are located inside the clinic.

Apart from these, consumables required for obstetrics and examination are also placed in the clinic.

There is also a hot and cold water supply system, water treatment system, waste water tank, hydrophore, generator. Generator booster, water tanks are located at the bottom of the vehicle.

There is a 2 meter wide awning/canopy on one side of the vehicle. The entrance of the clinic is made from the back and the automatic door with shutter system is made as an outer door. Apart from this, the inner door is the automatic operating room door with a card/password.

There is a hanger closet at the entrance and a computer desk with sockets next to it. Right next to it, there are counters and cabinets with drawers in the clinic.

For more detailed information & special design requests contact us via the contact information.



Positioned strategically within the interior of the Container Maternity Unit is a drug preparation station, meticulously crafted from PVC antibacterial material. This bench is thoughtfully situated near the central partition of the unit, ensuring convenient access. Additionally, an overhead cabinet system graces the bench, facilitating secure storage of medicines and needles.

The flooring of the maternity examination unit is adorned with an epoxy coating, rendering it effortlessly cleanable and hygienic.

The ceiling is adorned with 4 LED lights and 2 spot lamps, creating an illuminated and soothing atmosphere.

Within the unit, two 220v output sockets cater to devices requiring such power, with a total capacity of 2000w.

Should the need arise, the incorporation of a mobile ultrasound device can be seamlessly integrated into the unit.

Furthermore, external 220v power can be harnessed to energize the unit during stationary periods.

In response to specific requests, the option to include a discreet, low-emission 3kw gasoline generator adjacent to the right door is available.

Essential medical tools including a sphygmomanometer and birth control kits find their dedicated space within the unit. Moreover, upon request, blood test equipment can be furnished to augment the unit's capabilities.





# OUR WORKS





## **TECHNICAL EQUIPMENTS**

- 1 Mobile Health Care Vehicle
- 2 Front beacon lamp with LED
- 3 Back beacon lamp with
- 4 Siren Announce Device
- 5 Patient's Cabin Lighting (LED)
- 6 Patient's Cabin Lighting (Spot LED)
- 7 Air Condition
- 8 Electrical System
- 9 Digital Control Panel
- 10 Rectifier
- 11 Rectifier Cable
- 12 Socket
- 13 Inverter 1000W
- 14 Outlet 220V
- 15 Fuse Circuit Breaker
- 16 Led Lighting



The technical equipment in our ambulances is designed to meet the highest standards of quality and functionality. Please note that these equipment are subject to regular updates and advancements in line with evolving medical technologies. We are committed to ensuring that our ambulances are equipped with the most up-to-date and state-of-the-art technology to provide the best possible care.

## **MEDICAL EQUIPMENTS**

- 1 Delivery Bed
- 2 Sterile Supplies
- 3 Medical Lighting
- 4 Monitoring Equipment
- 5 Intravenous (IV) Equipment
- 6 Infant Resuscitation Equipment
- 7 Surgical Tools
- 8 Emergency Medications
- 9 Incubator/Warming Equipment
- 10 Cord Clamps and Scissors
- 11 Delivery Kits
- 12 Suction Devices
- 13 Maternal and Newborn Care Supplies
- 14 Waste Disposal
- 15 Emergency Transport Equipment
- 16 Privacy Screens
- 17 Hand Hygiene and Sanitization Supplies
- 18 Emergency Obstetric Kit



**Ambulancemed**