

AMBULANCEMED

AmbulanceMed, who combines technical expertise with the commercial ones and maintains with national/international experience, has led the delivery of thousands of ambulances until today.

FACTORY - TURKEY

info@ambulancemed.com
sales@ambulancemed.com

+90 554 969 50 08

DAĞYAKA MAH. 2018 SK. NO:23
AKSAN TİCARET MERKEZİ
KAHRAMANKAZAN / ANKARA

FACTORY - DUBAI / UAE

info@ambulancemed.com
dubai@ambulancemed.com

+971 56 332 6520

X3WR+Q57, Mina Jebel Ali, N506
Jebel Ali Free Zone, Dubai, UAE

OFFICE - DUBAI / UAE

info@ambulancemed.com
dubai@ambulancemed.com

+90 554 969 50 08

Unit no. 604-605, 6F Opal Tower,
Business Bay/ Dubai/ 128846
Registration No: ICC20210433

OFFICE / INTERNATIONAL

info@ambulancemed.com
export@ambulancemed.com

+90 554 969 50 08

ATB İş Merkezi, G Blok No: 172
Yenimahalle / Ankara



Ambulancemed

FIAT DUCATO GYNECOLOGY VEHICLE

Rapid Response. Unmatched Care
Ducato Gynecology Vehicle.

 **MOBILE GYNOCLOGY CLINIC**





WELCOME **INTRODUCTION**

Step into the world of AmbulanceMed and experience the epitome of excellence in emergency medical services. Our visionary company proudly presents the exquisitely designed Ducato Gynecology Vehicle, meticulously crafted to redefine the standards of ambulance technology. In this captivating brochure, we invite you to embark on a detailed exploration of the AmbulanceMed masterpiece, where form seamlessly blends with function. Discover the unmatched sophistication of our cutting-edge features, engineered to deliver exceptional patient care and ensure swift, efficient emergency response. Join us as we push the boundaries of innovation and revolutionize the future of ambulance solutions.

TABLE OF CONTENT

You can explore a wealth of information about Ducato Gynecology Vehicle within this brochure. Please refer to the following page numbers to find detailed content on each topic

01. [Fiat Ducato Gynecology Vehicle](#)

02. [Vehicle Interior](#)

03. [Table Of Content Goes Here](#)

04. [Technical Specifications](#)

05. [Medical Specifications](#)

06. [Gallery](#)

07. [Technical Equipments](#)

08. [Medical Equipments](#)

You are about to embark on an informative journey delving deep into the realm of Gynecology Vehicle

FIAT DUCATO

GYNECOLOGY VEHICLE

At Ambulancemed, we take great pride in our commitment to transforming Fiat Ducato brand cars into top-of-the-line gynecology vehicle, and our latest creation, the Ducato Gynecology Vehicle, is no exception. Manufactured at our state-of-the-art facilities in Dubai and Ankara, this vehicle is designed to provide unparalleled medical care, ensuring optimal sanitary conditions and comfort for both patients and medical professionals.

The essence of the Gynecology Vehicle's design lies in its patient compartment, meticulously crafted to provide the highest standards of medical care. Understanding the critical nature of gynecological services, we have ensured ample space for comfortable patient treatment. Every inch of the patient compartment is thoughtfully designed, allowing for efficient storage of specialized medical supplies, making them readily accessible during critical gynecological emergencies.

Our vehicle. Medical devices, vehicle connection mounting elements is manufactured as being appropriate with EN 1789+A2 standard.



As a leading manufacturer of medical vehicles, we are dedicated to serving not only local communities but also extending our reach worldwide. Our Ducato Gynecology Vehicle is exported to various countries, ensuring that people around the globe have access to top-tier gynecological medical services.

The Ducato Gynecology Vehicle exemplifies our unwavering commitment to excellence, safety, and innovation. With its unmatched patient care facilities, ergonomic design, efficient lighting and warning systems, and revolutionary WAS assistance systems, it stands at the forefront of the medical vehicle industry. Together, we strive to elevate gynecological medical care and make a significant impact on countless lives across the globe. Choose Ambulancemed - your trusted partner in delivering exceptional medical vehicle solutions.

■ Efficient Illumination and Special Warning System

The Ducato Gynecology Vehicle is equipped with an energy-efficient lighting system, providing ample illumination for gynecological procedures without excessive power consumption. Our special warning system is meticulously designed to garner attention while being conscious of energy use. The combination of effective lighting and warning systems ensures a safe passage through traffic, allowing the vehicle to reach its destination swiftly during critical gynecological situations.

At Ambulancemed, we believe in equipping our medical vehicles with the latest technologies. The Ducato Gynecology Vehicle comes equipped with the cutting-edge WAS assistance systems, including Door Assist, Multi-Load Assist, SanSafe, and Smart Check Mobile. These state-of-the-art technologies simplify the daily work of the medical team, streamlining operations and enhancing overall efficiency during gynecological emergencies and procedures.

- The Fiat Ducato Gynecology Vehicle by Ambulancemed is not only a testament to our dedication to safety and comfort but also a platform to showcase the latest technological medical devices, revolutionizing the landscape of gynecological healthcare.



Ambulancemed





VEHICLE INTERIOR

At AmbulanceMed, we take pride in equipping the Fiat Ducato Gynecology Vehicle interior with full technological advancements, keeping pace with the demands of the modern era. Our company is dedicated to providing the highest standard of gynecological healthcare, and we have meticulously designed the vehicle to meet those standards.

The interior is thoughtfully furnished with cutting-edge and tailored equipment, ensuring that it is fully equipped to handle any gynecological situation. From advanced medical devices specific to gynecology to user-friendly layouts, every detail has been carefully considered to meet the specific needs of gynecological medical professionals and patients alike. With AmbulanceMed, you can trust that our technologically advanced Fiat Ducato Gynecology Vehicle will deliver optimal care for women's health when it matters



Ambulancemed



GENERAL DESCRIPTIONS

FIAT DUCATO GYNECOLOGY VEHICLE

Mobile Gynecology Clinic could be mounted on a van or trailer. It was designed in order to make gynecological examine or make outpatient treatment to patients who live in rural areas and have limited access to the hospital.

The Mobile Gynecology Clinic Consists of Three Parts. There is an input area looked from the backside of the vehicle that has a space covered with a thick curtain inside it in order to enable patients to dress up and undress. A gynecological examination table, which is a full-function, is placed smoothly in the middle of the vehicle. In fact, that table allows making urologic and minor surgical operations.

- The backrests of the table can be arranged in terms of its height by an electric motor.
- A hand control located in the patient cabinet for the gynecological table.
- It has a shock absorber mechanism
- Can be taken the position of Trendelenburg
- Covered with a particular substance
- Durable
- Easily cleaned
- Consist of back, foot and seat part
- All the components are modular and can be disassembled at any time
- Contains foot and arm support elements in order to facilitate the examination
- A waste cuvette is based under the desk
- A doctors' seat, which has rotary type ground, is located on the left side of the desk



It is placed that a drug preparation bench, made from PVC antibacterial substance, near the intermediate partition of the vehicle. Moreover, there is a cabinet system on the upper side of the bench in order to store medicines and needles.

The floor of the gynecological examination ambulance was covered with an epoxy material that can be cleaned easily.

The ceiling is lighted with 4 LEDs and 2 spot lamps.

There are two 220v output sockets in the vehicle for the devices working with 220v. It enables the devices that 2000w in total to work.

It is mounted that the mobile ultrasound device into the vehicle on the condition that asked for.

It is feasible to provide electricity to the vehicle by using 220v cable that came from outside while the vehicle is not in motion.

A 3kw silent gasoline type generator, which has a minimum level of smoke emission, can be mounted on the next to the right door in case asked for.

There are sphygmomanometer and birth control sets placed in the vehicle; in addition, blood test equipment can be provided if requested.





OUR WORKS



TECHNICAL EQUIPMENTS

Mobile Health Care Vehicle Technical Details	ORIGIN	STANDARD	QTY
Front beacon lamp with LED	Turkey	ECR-65	1
Back beacon lamp with LED + projector	Turkey	ECR-65	1
Siren Announce Device (5 stage)	Turkey	ECR-65	1
Side warning lights with LED 12 volt	Turkey	ECR-65	6
Patient's Cabin Lighting (LED)	Turkey	CE	2
Patient's Cabin Lighting (Spot LED)	Holland	CE	2
Air Condition	Turkey	CE / ISO 9001	1
Electrical System	Turkey	CE / ISO 9001	1
Digital Control Panel	Turkey	CE / ISO 9001	1
Rectifier	Turkey	CE / ISO 9001	1

TECHNICAL EQUIPMENTS

Rectifier Cable	Turkey	CE / ISO 9001	1
Socket	Turkey	CE / ISO 9001	1
Inverter 1000W	Turkey	CE / ISO 9001	1
Outlet 220V	Turkey	CE / ISO 9001	1
Fuse Circuit Breaker	Turkey	CE / ISO 9001	1
Led Lighting	Turkey	CE / ISO 9001	1
Stretcher Spotlight	Turkey	CE / ISO 9001	1
Electronic Siren and Light Bar	Turkey	CE / ISO 9001	1
Back Top Lightbar	Turkey	CE / ISO 9001	1



MEDICAL EQUIPMENTS

MEDICAL EQUIPMENTS

Mobile Health Care Vehicle Medical Equipments Details	ORIGIN	STANDARD	QTY
Mobile Health Care Vehicle Main Stretcher	Turkey	CE / EN 1865 / EN 1789 / ISO 9001	1
Mobile Health Care Vehicle Chair Stretcher	Turkey	CE / EN 1865 / EN 1789 / ISO 9001	1
Mobile Health Care Vehicle Scoop Stretcher	Turkey	CE / EN 1865 / EN 1789 / ISO 9001	1
Mobile Health Care Vehicle Spine board Stretcher	Turkey	CE / EN 1865 / EN 1789 / ISO 9001	1
Mobile Health Care Vehicle Air Split Set	German/Turkey	CE / EN 1865 / EN 1789 / ISO 9001	1
Mobile Health Care Vehicle Cervical Coller	German	CE / EN 1865 / EN 1789 / ISO 9001	1
Mobile Health Care Vehicle Birth Set	German	CE / EN 1865 / EN 1789 / ISO 9001	1
Mobile Health Care Vehicle Burn Set	German	CE / EN 1865 / EN 1789 / ISO 9001	1
Mobile Health Care Vehicle Portable Suction Set	German/Turkey	CE / EN 1865 / EN 1789 / ISO 9001	1
Mobile Health Care Vehicle Defibrillator AED	German	CE / EN 1865 / EN 1789 / ISO 9001	1
Mobile Health Care Vehicle Transportation Ventilator	German	CE / EN 1865 / EN 1789 / ISO 9001	1

Mobile Health Care Vehicle Reanimation Set	German/Turkey	CE / EN 1865 / EN 1789 / ISO 9001	1
Mobile Health Care Vehicle Fixed Spyhg	German	CE / EN 1865 / EN 1789 / ISO 9001	1
Mobile Health Care Vehicle Portable Spyhg	German	CE / EN 1865 / EN 1789 / ISO 9001	1
Mobile Health Care Vehicle Stethoscope	German	CE / EN 1865 / EN 1789 / ISO 9001	1
Mobile Health Care Vehicle Corpse Bag	Turkey	CE / EN 1865 / EN 1789 / ISO 9001	1
Mobile Health Care Vehicle Basic Medical Equipment Bag	German/Turkey	CE / EN 1865 / EN 1789 / ISO 9001	1
Mobile Health Care Vehicle Oxygen and Vacuum Systems	Turkey	CE / EN 1865 / EN 1789 / ISO 9001	1
Mobile Health Care Vehicle Oxygen Cylinder 10Ltx2	Turkey	CE / EN 1865 / EN 1789 / ISO 9001	1
Mobile Health Care Vehicle Oxygen Outlet	Turkey	CE / EN 1865 / EN 1789 / ISO 9001	1
Mobile Health Care Vehicle Oxygen Regulator & Flowmeter	Turkey	CE / EN 1865 / EN 1789 / ISO 9001	1



Ambulancemed